



SeeStar[®]

Guiding Device



**Optimize efficiency during your
procedures with SeeStar[®]**

Precise
Accurate
Reliable

AprioMed

SeeStar® Guiding Device

- Visualize and adjust the planned CT-guided biopsy before entering the needle.
- Reduce the number of needle manipulations¹.
- Support the needle during superficial punctures.
- Maintain the direction of the needle during puncture.



SeeStar® aids the physician by producing an artifact in the CT image, indicating the intended needle path.

AprioMed provides tools that enable physicians to optimize precision and control during the most challenging procedures. For more information about AprioMed's Interventional Radiology products please visit www.apriomed.com.

Bonopty® Bone Biopsy System



SeeGrid® Position Indicator



Before using SeeStar® read the Instructions for Use which accompany the product. SeeStar® is a registered trademark of AprioMed AB, Uppsala, Sweden. SeeStar® is manufactured under one or more of the following patents: EP 1534153, AU 2003256201, SE 524 042. Other patents pending.
¹ "CT Fluoroscopy Guided Punctures Performed with a Freehand Technique and in Combination with a New Needle Guiding Device" Anders Magnusson, Jon Isacson, Pär Dahlman, Allina Dimopoulou, Lars Jangland, Maria Lönnemark, Eva Radecka & Hans Raland, Depts. of Radiology and Medical Physics, University Hospital, Uppsala, Sweden, presented at ESTI 2007, Athens Greece.

ORDERING INFORMATION

Product number	Product Description	Minimum Order Quantity
13314	SeeStar® 14 gauge (green) D, 2,2 mm 0,087 inch, suitable for needle sizes 14 and 15 gauge	10 units
13317	SeeStar® 17 gauge (pink) D, 1,6 mm 0,063 inch, suitable for needle sizes 17 and 18 gauge	10 units
13320	SeeStar® 20 gauge (yellow) D, 1,0 mm 0,039 inch, suitable for needle sizes 20 and 21 gauge	10 units

29003 Rev. 04 2009-05

Head Office:
AprioMed AB
Virdings Allé 28
SE-754 50 Uppsala, Sweden
Phone: +46 (0)18 430 14 40
Free Fax (Int'l): +46 (0)770 339972

E-mail: info@apriomed.com
www.apriomed.com

US Sales Office:
AprioMed, Inc.
2 Palmer Dr, Suite 1a
Londonderry, NH 03053, USA
Phone: +1 603 421 0875
Fax (US toll free): 866 239 1797