

SeeGrid[®]

Position Indicator



**Optimize efficiency during your
procedures with SeeGrid[®]**

Precise
Accurate
Reliable

AprioMed

SeeGrid® Position Indicator

- A marker for all procedures that require position information to be transferred from the CT image to the patient's skin.
- Increases position transfer accuracy.
- Good visibility in all x-ray based image modalities.



SeeGrid® is made from soft flexible plastic filled with barium to ensure good visibility in all x-ray based image modalities.

AprioMed provides tools that enable physicians to optimize precision and control during the most challenging procedures. For more information about AprioMed's Interventional Radiology products please visit www.apriomed.com.

Bonopty® Bone Biopsy System



SeeStar® Guiding Device



Before using SeeGrid® read the Instructions for Use which accompany the product. SeeGrid® is a registered trademark of AprioMed AB, Uppsala, Sweden. SeeGrid® is manufactured under one or more of the following patents: SE 526 587, US 7,086,172, AU 2004202432. Other patents pending.

ORDERING INFORMATION

| Product number | Product Description | Minimum Order Quantity |
|----------------|--------------------------------------|------------------------|
| 13200 | SeeGrid® Position Indicator, 21x8 cm | 10 units |

29095 Rev. 06 2010-11

Head Office:
AprioMed AB
Virdings Allé 28
SE-754 50 Uppsala, Sweden
Phone: +46 18 430 14 40
Free Fax (Int'l): +46 770 339972

E-mail: info@apriomed.com
www.apriomed.com

US Sales Office:
AprioMed, Inc.
2 Palmer Dr, Suite 1a
Londonderry, NH 03053, USA
Phone: +1 603 421 0875
Fax (US toll free): 866 239 1797