



AprioMed

DuoGuide[®] II
Guidewire

TaperTip[™]
Dilator



Optimize tactility and efficiency during the first stage of the nephrostomy procedure with DuoGuide[®] II and TaperTip[™]

Precise
Accurate
Reliable

DuoGuide® II Guidewire

TaperTip™ Dilator


Designed to make the first stage of the nephrostomy procedure simpler and less painful

DuoGuide® II Guidewire

- Catheter support without compromising formability
- Achieve accurate placement through visualization
- Introduce through a 0.9 mm needle


TaperTip™ Dilator

- Achieve dilation in a single step
- Lessen pain with the long tapered tip
- Reduce risk of penetrating the collecting system



AprioMed provides tools that enable physicians to optimize precision and control during the most challenging procedures. For more information about AprioMed's Interventional Radiology products please visit www.apriomed.com

Bonopty® Bone Biopsy System




SeeStar® Guiding Device



SeeGrid® Position Indicator



 Prior to using DuoGuide® II or TaperTip™, read the Instructions for Use accompanying the products.

DuoGuide® II and TaperTip™ are not available for sale in the United States.
DuoGuide® II and TaperTip™ are trademarks of AprioMed AB, Uppsala, Sweden.

ORDERING INFORMATION

Product number	Product description	Minimum order quantity
13002 NEW	DuoGuide® II Guidewire Tip Flexibility Std super-elastic with shapeable tip <ul style="list-style-type: none"> • Ø 0.46 mm / 0.018" • 900 mm length • 45 mm soft platinum tip • Nitinol shaft 	5 units
13004 NEW	DuoGuide® II Guidewire Tip Flexibility High super-elastic with shapeable tip <ul style="list-style-type: none"> • Ø 0.46 mm / 0.018" • 900 mm length • 80 mm soft platinum tip • Nitinol shaft 	5 units
13100	TaperTip™ Dilator radioopaque, compatible with 0.018" guidewire <ul style="list-style-type: none"> • Inner Ø 0.51 mm • Outer Ø 3.3 mm • 30 mm tapered tip • 240 mm length 	10 units

29158 Rev. 01 2011-03

Head Office:
AprioMed AB
Virdings Allé 28
SE-754 50 Uppsala, Sweden
Phone: +46 18 430 14 40
Free Fax (Int'l): +46 770 339972

E-mail: info@apriomed.com
www.apriomed.com

AprioMed